


NOAH-NMR Supersequences with Nested Acquisition for Small Molecules

Product used : Nuclear Magnetic Resonance (NMR)


NOAH (**N**MR by **O**rdered **A**cquisition using ¹**H**-detection) [1] is a group of nested NMR experiments combining several conventional two-dimensional (2D) NMR pulse sequences, such as COSY, HSQC and HMBC, into one supersequence. Therefore, two or more 2D NMR data can be obtained from a single NOAH experiment. By using a single relaxation delay, the NOAH method significantly reduces the total data collection time and increases the throughput of an NMR instrument in structure elucidation of small organic molecules.


NOAH Supersequence

Figure 1 is a schematic representation of a conventional 2D NMR measurement, and Figure 2 shows the principle of a NOAH supersequence. All conventional 2D NMR sequences contain a recovery delay (relaxation delay) in the order of seconds, which is the longest part of the sequence. If two or more conventional 2D experiments are run in a sequence, each experiment needs to have a recovery delay (Figure 1). In the NOAH supersequence, only a single recovery delay is employed at the very beginning of the nested sequence (Figure 2). As a result, the total data collection time of the NOAH experiment is much shorter in comparison to the classical 2D data collection scheme.










[1] Ē. Kupče, T. D. W. Claridge, Angew. Chem. Int. Ed. 2017, 56, 11779.

JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

The Example of NOAH-3 Experiment

In Figure 3 (a), there is a NOAH-3 supersequence, that combines HMBC, multiplicity-edited HSQC and COSY. The ¹H-¹³C HMBC (b), ¹H-¹³C multiplicity-edited HSQC (c) and ¹H-¹H COSY (d) spectra were collected in 10 minutes. For comparison, it would take nearly 30 minutes to collect the spectra conventionally.



Copyright © 2019 JEOL Ltd.

Local office

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

3-1-2 Musashino Akishima Tokyo 196-8558 Japan Sales Division Tel. +81-3-6262-3560 Fax. +81-3-6262-3577 www.jeol.com ISO 9001 · ISO 14001 Certified

OREA • MALAYSIA • MEXICO • RUSSIA • SCANDINAVIA

• AUSTRALIA & NEW ZEALAND • BELGIUM • BRAZIL • CANADA • CHINA • EGYPT • FRANCE • GERMANY • GREAT BRITAIN & IRELAND • ITALY • KOREA • MALAYSIA • MEXICO • RUSSIA • SCANDINAVIA • SINGAPORE • TAIWAN • THE NETHERLANDS • USA

JEOL Ltd.

JEOL