

Qualitative Analysis of Calcium Oxalate Using TG-MS System

Product used: Mass Spectrometer (MS)

Overview

A thermogravimetry/differential thermal analysis (TG/DTA) system is used to measure weight changes and relative temperatures of samples under programmed heat conditions. A system combining TG/DTA with mass spectrometry (MS) is called a TG-MS system, and is considered one of the best qualitative and quantitative analysis systems for inorganic materials because it can measure sample conditions, evolved gas species, and gas volume simultaneously. In this application note, we introduce the basic application of analyzing calcium oxalate using the TG-MS “STA2500 Regulus” system (NETZSCH) and the GC/QMS “JMS-Q1500GC” system (JEOL).

Equipment and software

The TG/DTA and gas chromatography/MS (GC/MS) are connected by an inactivated fused-silica capillary tube in the TG transfer line, which is kept at a high temperature. The evolved gases from TG/DTA are introduced into the MS through this tube, which allows the system to measure both MS and TG signals simultaneously. Joint-analysis TG/DTA–GC/MS software can import TG/DTA and MS data, respectively. Consequently, it is possible to correlate weight change results with qualitative-analysis results for evolved gases.





Figure 1. Device and PC software connections



Experiment

Calcium oxalate hydrate samples of 5 mg and 20 mg were weighed with an electronic balance. Table 1 shows the TG-MS measurement conditions.

Table 1. TG-MS measurement conditions.

TG conditions		MS conditions	
Furnace temp.	60°C → 20°C/min → 1000°C	Ion source temp.	250°C
Transfer line temp.	300°C	Interface temp.	300°C
Atmospheric gas flow	He, 100 mL/min	Ionization mode	EI, 70 eV
Split ratio	100:1	Ionization current	50 μA
GC conditions		Relative EM voltage	+200 V
Oven temp.	300°C	Measurement mode	Scan
Column	Inactivated fused silica capillary tube (metal), 5 m × 0.25 mm i.d	Scan range	m/z 10–100

Results

Figure 2 shows measurement results for the 20 mg sample. TG curve (= weight change), DTA curve, and total ion current chromatograms (TICCs) are respectively shown as green, blue, and red lines. On the TG/DTA curve, the weight loss from the endothermic reaction was detected at 190°C, 500°C, and 780°C. The TICC shows three evolved gas peaks at the same temperature as the weight losses. Figure 3 shows the mass spectra of the TICC peaks [1]–[3] from Figure 2. Peaks [1], [2], and [3] were assigned as H₂O, CO/CO₂, and CO₂, respectively, by NIST Library search.




Figure 2. TG/DTA curves and TICC for calcium oxalate.





Figure 3. Mass spectra of TICC peaks.

The following chemical-reaction formula shows the thermal-decomposition process of calcium oxalate. The TG-MS measurement results suggest the following chemical reaction occurred:



Repeatability

Table 2 shows the extracted-ion chromatogram (EIC) peak area (per weight) of calcium oxalate 5 mg (n = 5) and its relative standard deviation (RSD). The RSD values of EIC peaks for the three evolved gases are within 3%.

Table 2. Area of EIC peak.

	[1]H ₂ O (m/z 18)	[2]CO (m/z 28)	[3]CO ₂ (m/z 44)
1	72522896	140532864	384935646
2	74839998	144486163	393943250
3	73120172	142488903	385202014
4	73118409	141937036	385451041
5	69352941	135562682	367173088
RSD	2.8%	2.4%	2.6%

Copyright © 2023 JEOL Ltd.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.



JEOL Ltd.

3-1-2 Musashino Akishima Tokyo 196-8558 Japan Sales Division Tel. +81-3-6262-3560 Fax. +81-3-6262-3577
www.jeol.com ISO 9001 • ISO 14001 Certified

• AUSTRALIA & NEW ZEALAND • BELGIUM • BRAZIL • CANADA • CHINA • EGYPT • FRANCE • GERMANY • GREAT BRITAIN & IRELAND • ITALY • KOREA • MALAYSIA • MEXICO • RUSSIA • SCANDINAVIA
• SINGAPORE • TAIWAN • THE NETHERLANDS • USA

Local office

